

Table 1A. Criteria for the classification of pathogenic mt-tRNA variants. The ACMG criteria are listed on the left column and the new criteria for mt-tRNA variants are on the right.

Evidence	ACMG criteria	Criteria for mt-tRNA variants
Very strong	<b>PVS1</b> Null variant in a gene where loss of function is a known mechanism of disease.	<b>PVS1</b> Not applicable.
Strong	<p><b>PS1</b> Same amino acid change as a previously established pathogenic variant regardless of nucleotide change.</p> <p><b>PS2</b> <i>De novo</i> in a patient with disease and no family history.</p> <p><b>PS3</b> Well-established <i>in vitro</i> or <i>in vivo</i> functional studies supportive of a damaging effect on the gene or gene product.</p> <p><b>PS4</b> The prevalence of the variant in affected individuals is significantly increased compared to the prevalence in controls.</p>	<p><b>PS1</b> Not applicable.</p> <p><b>PS2</b> Present at <math>\geq 5\%</math> heteroplasmy, in <math>\geq 2</math> different tissues of the affected individual but 0% in asymptomatic mother. If mother's sample is unavailable, and 0% in other asymptomatic matrilineal relatives (such as proband's siblings, proband's maternal grandmother), it will be downgraded to PM9. Note: The percentage of heteroplasmy should be analyzed by dependable, clinically validated method, such as deep NGS using one-piece Long Range PCR product of the circular mtDNA template. The sensitivity of the methodology must be able to distinguish true zero from background zero.</p> <p><b>PS3</b> Well-established <i>in vitro</i> or <i>in vivo</i> functional studies supportive of a damaging effect on the mitochondrial function. These could be transmitochondrial cybrid assays, ETC, OCR, ATP synthesis, mtDNA copy number, COX deficient fibers, single fibers, etc, at a diagnostic levels (depending on assays and tissue types), and correlate with percentage of heteroplasmy. If not correlate with percentage of heteroplasmy, it will be downgraded to PM10.</p> <p><b>PS4</b> The prevalence of the variant in affected individuals is significantly increased compared to the prevalence in controls.</p>

		<p><b>PS5</b> Rare variants previously reported as pathogenic.</p> <p>Note: Not all reports are reliable, especially those old ones published before the application of deep NGS to sequence mtDNA. Please review the literature cautiously before use this criterion.</p>
<b>Moderate</b>	<p><b>PM1</b> Located in a mutational hot spot and/or critical and well-established functional domain without benign variants</p> <p><b>PM2</b> Absent from controls in Exome Sequencing Project, 1000 Genomes or ExAC.</p> <p><b>PM3</b> For recessive disorders, detected in <i>trans</i> with a pathogenic variant.</p> <p><b>PM4</b> Protein length changes due to in-frame deletion/insertions in a non-repeat region or stop-loss variants.</p> <p><b>PM5</b> Novel missense change at amino acid residue where a different missense change determined to be pathogenic has been seen before.</p> <p><b>PM6</b> Assumed <i>de novo</i>, but without confirmation of paternity and maternity.</p>	<p><b>PM1</b> Variants that cause anticodon swap.</p> <p><b>PM2</b> Absent from database, e.g., mtDB and MitoMap, and absent or low heteroplasmy (&lt;5%) in asymptomatic mother, compared to that in the proband. If mother's sample is unavailable, it will be downgraded to PP7.</p> <p>Note: The percentage of heteroplasmy should be analyzed by dependable, clinically validated method, such as deep NGS, using one-piece Long Range PCR product of the circular mtDNA template.</p> <p><b>PM3</b> Not applicable.</p> <p><b>PM4</b> Not applicable.</p> <p><b>PM5</b> Same position nucleotide change of a previously well-established pathogenic variant to a different nucleotide. Example: m.3243A&gt;T vs m.3243A&gt;G.</p> <p><b>PM6</b> Not applicable. See PS2</p> <p><b>PM7</b> MitoTIP prediction score &gt;16.0.</p>

		<p><b>PM8</b> Heteroplasmy (<math>\geq 5\%</math>) among different tissues of an affected individual correlates with clinical or biochemical phenotypes. Example: The heteroplasmy of a patient's muscle is at 10%, while the heteroplasmy in blood is at 3%.</p> <p><b>PM9</b> At least two independent families, or two matrilineal family members from one family demonstrate correlation of heteroplasmy (<math>\geq 5\%</math>) with clinical or biochemical phenotypes.</p> <p><b>PM10</b> Well-established <i>in vitro</i> or <i>in vivo</i> functional studies supportive of a damaging effect on the mitochondrial function. These could be transmitochondrial cybrid assays, ETC, OCR, ATP synthesis, mtDNA copy number, COX deficient fibers, ragged red fibers, etc, at a diagnostic levels (depending on assays and tissue types), <b>but not correlate with percentage of heteroplasmy.</b></p>
<b>Supporting</b>	<p><b>PP1</b> Co-segregation with disease in multiple affected family members in a gene definitively known to cause the disease.</p> <p><b>PP2</b> Missense variant in a gene that has a low rate of benign missense variation and where missense variants are a common mechanism of disease.</p> <p><b>PP3</b> Multiple lines of computational evidence support a deleterious effect on the gene or gene product.</p> <p><b>PP4</b> Patient's phenotype or family history is highly specific for a disease with a single genetic etiology.</p> <p><b>PP5</b> Reputable source recently reports variant as pathogenic but the evidence is not available to the laboratory to perform an independent evaluation.</p>	<p><b>PP1</b> Not applicable. Upgrade to PM8 or PM9 when co-segregate not only with disease but also with <math>&gt;5\%</math> heteroplasmy.</p> <p><b>PP2</b> Not applicable.</p> <p><b>PP3</b> The range of MitoTIP prediction score is within [12.5-16].</p> <p><b>PP4</b> Patient's phenotype or family history is highly specific for a mitochondrial disease with a single genetic etiology.</p> <p><b>PP5</b> Not applicable.</p> <p><b>PP6</b> Heteroplasmy in an affected proband <math>\geq 5\%</math>.</p>

---

**PP7** Absent from database, e.g., mtDB and MitoMap, and is heteroplasmic.

---

ETC: electron transport chain, OCR: oxygen consumption rate, COX: cytochrome C oxidase. mtDB: human mitochondrial genome database.

Table 1B. Criteria for the classification of benign variants.

Evidence	ACMG criteria	Criteria for mt-mRNA variants
<b>Stand-alone</b>	<b>BA1</b> Allele frequency is above 5% in Exome Sequencing Project, 1000 Genomes, or ExAC.	<b>BA1</b> Top-level haplogroup defining variants.
<b>Strong</b>	<p><b>BS1</b> Allele frequency is greater than expected for disorder.</p> <p><b>BS2</b> Observed in healthy adult individual for a recessive (homozygous), dominant (heterozygous), or X-linked (hemizygous) disorder with full penetrance expected at an early age.</p> <p><b>BS3</b> Well-established <i>in vitro</i> or <i>in vivo</i> functional studies shows no damaging effect on protein function or splicing.</p> <p><b>BS4</b> Lack of segregation in affected members of a family.</p>	<p><b>BS1</b> Reported in public databases (e.g., MitoMap or mtDB) or literatures as polymorphism.</p> <p><b>BS2</b> Found homoplasmic in more than three unrelated healthy adults.</p> <p><b>BS3</b> Not applicable.</p> <p><b>BS4</b> Homoplasmy in both probands and at least 2 asymptomatic matrilineal family members.</p>
<b>Supporting</b>	<p><b>BP1</b> Missense variant in a gene for which primarily truncating variants are known to cause disease.</p> <p><b>BP2</b> Observed in <i>trans</i> with a pathogenic variant for a fully penetrant dominant gene/disorder; or observed in <i>cis</i> with a pathogenic variant in any inheritance pattern.</p> <p><b>BP3</b> In-frame deletions/insertions in a repetitive region without a known function.</p> <p><b>BP4</b> Multiple lines of computational evidence suggest no impact on gene or gene product.</p> <p><b>BP5</b> Variant found in a case with an alternate molecular basis for disease.</p>	<p><b>BP1</b> Not applicable.</p> <p><b>BP2</b> Not applicable.</p> <p><b>BP3</b> Not applicable.</p> <p><b>BP4</b> The MitoTIP prediction score is <math>\leq 10</math>.</p> <p><b>BP5</b> In the presence of a known pathogenic genetic cause unless there is evidence of more than one disease and clinically</p>

---

explained.

**BP6** Reputable source recently reports variant as benign but the evidence is not available to the laboratory to perform an independent evaluation.

**BP6** Found to be homoplasmic 1-3 times in private reputable laboratory databases.

**BP7** A synonymous (silent) variant for which splicing prediction algorithms predict no impact to the splice consensus sequence nor the creation of a new splice site AND the nucleotide is not highly conserved.

---

**BP7** Not applicable.